

Specification Sheet

Nitrile glove- powder free

Product Details

- XS, S, M, L, XL
- Physical measurement based on M size: L:9", ≈3mil, W:3.5g/pc(=/- 0.30 gram)
- blue
- smooth surface
- ambidextrous
- straight fingers
- micro-textured finger cots
- rolled cuff
- powder free

Material

- 100% nitrile

Norms & Guidelines

- CE according to Medical Device Directive 93/42/EEC Annex VII – Class I product
- PPE Category III - obligatory monitoring module C2
- Disposable product
- Latex-free
- Non sterile
- AQL 1.5
- ISO 13485:2016
- EN 455-1:2020, EN 455-2:2015, EN 455-3:2015
- UEReg425/2016, EN374-1:2016 + A1:2018, EN374-5:2016

Packaging

100 pcs/inner box, 1.000 pcs/carton

Code	Color	Size	Export per carton	Carton Dimension
80100	blue	XS, S, M, L, XL	1.000	28,4 x 22,8 x 22



Subject to modification without prior notice Oct., 2021

Dispo-Dent B.V. E-mail: info@dispo-dent.nl Tel: 0321-318140 Mobile: 0651418777