



·Hygiene ·Dental ·Medical

# Specification Sheet


## Coverall w/hood - White

### Product Details

- Size S - XXXL
- White
- Elasticated hood, cuff, waist and ankles
- One way zip
- Sewn seams

### Material

- Spunbond polypropylene
- Breathable, solid, light
- ≈ 35 gsm (±2 gsm)



### Norms & Guidelines

- CE according to Medical Device Directive 93/42/EEC + 2007/47/EC – Class I product
- Disposable product
- Latex-free
- Non sterile
- ISO 13485:2016

### Packaging

- 1 pc/polybag, 50 pcs/carton

### Size in cms

Code	Size	Color	Export per carton	Gross Weight per carton	Net Weight per Carton	Carton Dimension
10111	S	white	50	6,5 kg	5,8 kg	49 x 32 x 42
10112	M	white	50	6,5 kg	5,8 kg	49 x 32 x 42
10113	L	white	50	6,5 kg	5,8 kg	49 x 32 x 42
10114	XL	white	50	6,7 kg	6,0 kg	49 x 32 x 44
10115	XXL	white	50	6,7 kg	6,0 kg	49 x 32 x 44
10115X	XXXL	white	50	6,7 kg	6,0 kg	49 x 32 x 44

Subject to modification without prior notice

Dec., 2020

Dispo-Dent B.V. E-mail: [info@dispo-dent.nl](mailto:info@dispo-dent.nl) Tel: 0321-318140 Mobile: 0651418777