


**Specification Sheet**  
**PP isolation Gown - Blue**

**Product Details**

- size 130x140cm
- blue
- knitted cuffs
- ties in neck and waist



**Material**

- spunbond polypropylene
- solid, light
- breathable
- ≈ 23 gsm (±2 gsm)

**Norms & Guidelines**

- CE according to Medical Device Directive 93/42/EEC + 2007/47/EC – Class I product
- Disposable product
- Latex-free
- Non-sterile
- ISO 13485:2016

**Packaging**

10 pcs/polybag, 100 pcs/carton

Code	Color	Export per carton	Gross Weight per carton	Net Weight per Carton	Carton Dimension
30022	blue	100 pcs	6,7 kg	6,0 kg	52 x 28 x 33 cm



Subject to modification without prior notice

Jul., 2020