


Specification Sheet

PP isolation Gown - Yellow

Product Details

- size 130x140cm
- yellow
- knitted cuffs
- ties in neck and waist



Material

- spunbond polypropylene
- solid, light
- breathable
- ≈ 23 gsm (±2 gsm)

Norms & Guidelines

- CE according to Medical Device Directive 93/42/EEC + 2007/47/EC – Class I product
- Disposable product
- Latex-free
- Non-sterile
- ISO 13485:2016

Packaging

10 pcs/polybag, 100 pcs/carton

Code	Color	Export per carton	Gross Weight per carton	Net Weight per Carton	Carton Dimension
30024	yellow	100 pcs	6,7 kg	6,0 kg	52 x 28 x 33 cm



Subject to modification without prior notice

Jul., 2020