


**Specification Sheet**  
**PE coated isolation Gown - Yellow**

**Product Details**

- size 115x137cm
- yellow
- PE complete coating
- knitted cuffs
- ties in neck and waist



**Material**

- spunbond polypropylene+polyethylene
- solid, light
- moisture-resistant
- ≈ 30 gsm (±3 gsm)

**Norms & Guidelines**

- CE according to Medical Device Directive 93/42/EEC + 2007/47/EC – Class I product
- Disposable product
- Latex-free
- Non-sterile
- ISO 13485:2016

**Packaging**

10 pcs/polybag, 100 pcs/carton

Code	Color	Export per carton	Gross Weight per carton	Net Weight per Carton	Carton Dimension
30054	yellow	100 pcs	7,2 kg	6,5 kg	52 x 28 x 33 cm



Subject to modification without prior notice

Jul., 2020